PREVELLE[™] *silk*

MENTOR

3041 Skyway Circle North
Irving, TX 75038, USA

PREVELLE[™] Silk Dermal Filler

DESCRIPTION

PREVELLE Silk is a sterile, nonpyrogenic, viscoelastic, clear, colourless, transparent gel implant composed of cross-linked molecules of hyaluronan. Hyaluronan is a naturally occurring polysaccharide of the intercellular matrix in human tissues, including skin. Hyaluronan is chemically, physically and biologically identical in the tissues of all animal species.

COMPOSITION

Each 1 ml contains: hyaluronan gel 5,5 mg, lidocaine hydrochloride 3 mg, phosphate buffered saline, sodium chloride 9,3 mg, water for injection USP q.s.

MODE OF ACTION

PREVELLE Silk is injected into the dermal tissue to provide a space-occupying viscoelastic supplement for the intercellular matrix of the connective tissue. This viscosupplementation or augmentation of the dermal tissue can result in the temporary correction of skin contour deformities caused by wrinkles and depressed scars. Lidocaine in the formulation is intended to reduce pain during injections.

INDICATIONS AND USAGE

PREVELLE Silk is indicated for the temporary correction of cutaneous contour deformities (wrinkles, folds, scars) of the skin, particularly in cases of aging or degenerative lesions.

CONTRAINDICATIONS

PREVELLE Silk is contraindicated for patients with severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies. PREVELLE Silk is contraindicated for use in breast augmentation, or for implantation into bone, tendon, ligament, or muscle.

PREVELLE Silk contains trace amounts of gram positive bacterial proteins and is contraindicated for patients with a history of allergies to such material.

PREVELLE Silk contains lidocaine hydrochloride and must not be used in patients with known hypersensitivity to lidocaine or other local anaesthetics of the amide type.

WARNINGS

- PREVELLE Silk must not be injected into blood vessels. Introduction of PREVELLE Silk into the vasculature may occlude the vessels, and could cause infarction or embolisation.
- Use of PREVELLE Silk at specific sites in which an active inflammatory process (skin eruptions such as cysts, pimples, rashes, or hives) or infection is present, should be deferred until the underlying process has been controlled.

PRECAUTIONS

PREVELLE Silk is packaged for single patient use. Do not resterilise. Do not use if package is opened or damaged.

As with all percutaneous procedures, PREVELLE Silk use carries a risk of infection.

PREVELLE Silk must not be injected in patients with any acute or chronic skin disease in or near the intended area of correction.

The safety of PREVELLE Silk in patients with known susceptibility to keloid formation, hypertrophic scarring and pigmentation disorders has not been studied.

The safety of PREVELLE Silk use during pregnancy or in persons younger than 18 years of age has not been studied.

Patients who are using substances that can prolong bleeding, such as aspirin, non-steroidal anti-inflammatory drugs and warfarin may, as with any injection, experience increased bruising or bleeding at injection sites.

ANTICIPATED SIDE EFFECTS/ADVERSE EVENTS

The following types of adverse reactions may occur with any dermal medicinal products, such as PREVELLE Silk:

• Common reactions (> 1%):

- Patient reported local symptoms: Pain/tenderness and lumps/bumps (both 50% of patients); Swelling and

redness (both 48% of patients); Bruising (36% of patients); Itching (16% of patients); and Skin discoloured (9% of patients). In most cases the symptoms were mild in intensity and resolved within 3 days.

- Physician reported local adverse events: Injection site pain (4% of patients) which was mild or moderate in intensity and resolved in 2 days.

- Physician reported non-local adverse events: Gastroenteritis viral (2%), Onychomycosis (2%), Nasopharyngitis (2%) and skin irritation (2%).

• Uncommon reactions (> 0,1% but < 1%) and rare reactions (< 0,1%):

None were reported.

A controlled U.S. study showed use of PREVELLE Silk with lidocaine hydrochloride resulted in significantly less pain as compared to a hyaluronic acid based product without lidocaine hydrochloride.

INTERACTIONS

The safety of the combined use of PREVELLE Silk with other drugs and devices has not been established.

A medical history should be taken prior to treatment to ensure that there are no contraindications to the therapy. The patient's wrinkle and scar treatment areas should be properly identified and evaluated before treatment.

DIRECTIONS FOR USE

The clinical utility of PREVELLE Silk has been demonstrated in the correction of facial wrinkles, folds and scars (for example, acne and viral pox). Markedly indurated deficiencies and the superficial rhytides are difficult to correct. Depending on type of skin and lesion, best results are obtained in areas where these defects are readily distensible and where the correction can be visualised by manual manipulation (stretching) of the skin.

To assemble needle to syringe 1. Peel sealed cover off the needle guard. 2. Remove tip cap from syringe. 3. Attach needle to syringe and twist to secure. To assure proper needle attachment, use the needles provided or 30 gauge needles with similar needle guards. Fully seat hub of needle in syringe. Do not overtighten, as this may break the needle and/or dislodge the syringe. 4. Pull off the needle guard to expose needle.

Prior to treatment, patients should be advised about methods of administration, expected duration of correction, labeling precautions, possible treatment responses and adverse reactions.

Patients should also be advised that, in some cases, additional touch-up injections may be required to achieve appropriate correction.

The treatment area should be cleansed with alcohol or other antiseptic before the implant procedure. Prior to injecting PREVELLE Silk, depress the plunger rod until the

product flows out of the needle.

PREVELLE Silk is implanted through a narrow gauge needle (usually 30 G). The needle should be placed in the superficial to mid/deep dermis, depending upon the lesion being corrected. A linear threading technique, serial puncture injections, or a combination of the two have been used to achieve optimal results. Defects should be fully corrected at each treatment session. The defect should not be overcorrected. After PREVELLE Silk has been injected, the skin may be lightly massaged to conform to the contour of the surrounding tissue.

The rate and degree of subsidence of correction in the implanted area vary with patient, treatment site, and aetiology of the defect.

After the first treatment, one or more additional injections (at least one week apart) of PREVELLE Silk may be necessary to achieve the desired level of correction.

To sustain correction, repeat injections over time may be necessary (i.e., after 5 - 12 months). The need for reinjection may vary from site to site and is dependent upon a variety of factors, including anatomic location and the aetiology of the contour deficiencies and the technique of application.

Severely indurated defects, such as scars which initially resist distension, may require several treatment sessions before the desired correction is obtained. Such defects, however, tend to retain their correction over a longer period of time, so care should be taken not to overcorrect these sites.

HOW SUPPLIED

PREVELLE Silk is supplied in individual treatment syringes with two 30G 1/2 inch needles, and is packaged for single patient use, ready for injection (implantation). The syringe contents were sterilised by steam sterilisation and are non-pyrogenic. The volume markings on the syringe are for reference only.

1 x 0,75 ml

STORAGE DIRECTIONS

Store at room temperature, between 2°-30°C. Keep the product in the product box until used. DO NOT FREEZE. Protect against sunlight.

PRODUCT EVALUATION

Mentor requires that all complications from the use of this product be brought to the immediate attention of your local Mentor representative, who will be responsible for notifying Mentor. To contact Mentor directly, please telephone +31-71-5249600, fax +31-71-5217422 or e-mail to PEDepartment@mentormedical.nl.

**IF PACKAGE IS OPENED OR DAMAGED, DO NOT USE.
DO NOT USE IF SYRINGE CONTENTS ARE CLOUDY
OR DISCOLOURED**

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PREVELLE™ is a trademark and Mentor® is a registered trademark of Mentor Corporation.